

**MAX-LITE**<sup>®</sup>



## Ideal for long-term users and those with smaller ear canals

*MAX-LITE* (NRR 30) low pressure foam ear plugs offer long-term user comfort not found in other pre-shaped foam ear plugs. The non-irritating, non-allergenic self-adjusting foam recovers to fit virtually any ear canal.

Our revolutionary manufacturing process separates *MAX-LITE* from other ear plugs with a design that is extremely effective in reducing noise.

### **Comfortable**

There is virtually no pressure on the ear canal allowing for exceptional comfort.

### **Ideal for users with smaller ear canals**

The unique design of *MAX-LITE* ear plugs rolls extremely small.

### **Non-allergenic**

Soft, non-irritating, non-allergenic foam allows for wearer comfort.

### **NRR 30**

A highly rated disposable ear plug.

## Disposable Ear Plug fitting instructions

► Before fitting any ear plugs, make sure your hands are clean.

1. Hold the ear plug between your thumb and forefinger. Roll and compress the entire ear plug to a small, crease-free cylinder.



2. While still rolling, use your other hand to reach over your head and pull up and back on your outer ear. This straightens the ear canal, making way for a snug fit.

3. Insert the ear plug and hold for 20 to 30 seconds. This allows the ear plug to expand and fill your ear canal.



► **Test the fit.** In a noisy environment, and with earplugs inserted, cup both hands over your ears and release. You should not notice a significant difference in the noise level. If the noise seems to lessen when your hands are cupped over your ears, your ear plugs are probably not fitted properly. Remove and refit following instructions.



► Always remove ear plugs slowly, twisting them to break the seal. If you remove them too quickly, you could damage your ear drum.



► Foam ear plugs are disposable and not intended for re-use.



## Packaging options

MAX-LITE ear plugs are available packed in individual poly bags, and both corded and uncorded styles are available. The corded version of MAX-LITE offers greater convenience for users who frequently use and remove hearing protection. They're also more visible to help compliance supervisors ensure that employees are wearing hearing protection.



The Leight Ear Plug Dispenser is a permanent protection center that is easy to use, and buying ear plugs in bulk means you're not paying for unnecessary packaging.

## Attenuation Data

### MAX-LITE

ATTENUATION - USA (ANSI)									
TESTED ACCORDING TO ANSI S3.19-1974, R1979									
									NRR 30dB
Frequency (Hz)	125	250	500	1000	2000	3150	4000	6300	8000
Mean Attenuation dB(A)	33.5	33.6	36.0	37.5	39.4	42.5	43.9	43.7	45.2
Standard Deviation dB(A)	3.6	3.4	3.2	3.5	3.5	3.4	5.1	4.8	5.1

NRR 30 SNR 34CE SLC<sub>80</sub> Uncorded 27.0, Corded 27.5 CSA A(L)

## Ordering Information

- LPF-1 Without cord, packed in individual poly bags  
Case of 2,000 pair (10 dispenser boxes/200 pair per box)
- LPF-30 With poly cord, packed in individual poly bags  
Case of 1,000 pair (10 dispenser boxes/100 pair per box)
- LPF-30-P With white cotton cord, packed in biodegradable repulpable paper bags  
Case of 1,000 pair (10 dispenser boxes/100 pair per box)
- LPF-30-PB With poly cord, packed in biodegradable repulpable paper bags  
Case of 1,000 pair (10 dispenser boxes/100 pair per box)
- LPF-1-D Dispenser box  
Bulk packed for use in the Howard Leight dispenser ONLY  
Case of 2,000 pair (4 dispenser boxes/500 pair per box)

## HOWARD LEIGHT INDUSTRIES A Division of Bacou USA Safety, Inc.

7828 Waterville Road  
San Diego, CA 92173  
Tel (800) 327-1110  
Fax (800) 322-1330  
www.howardleight.com

